Synthesis Preloaded

PRELOADED ABERRATION CUTTING EDGE *sense of innovation*

FREE ASPHERIC LENS

Models SIP and SIPY (yellow)

EN INTRAOCULAR LENS

DESCRIPTION:
The Synthesis intraocular lens (IOL) is a single-piece ultraviolet (UV) absorbing posterior chamber IOL developed to replace the natural crystalline lens in pseudophakic adult patients. The Synthesis IOL has aspheric optics designed to be free of spherical aberration. The preloaded intraocular Synthesis lens is supplied with an Accuject Pro Delivery System to provide a sterile, controlled and touch-free method of delivery of the lens into the ocular bag. The Preloaded Synthesis IOL is supplied in 2 different models - Models SIP or SIPY. See label on the cardboard box for type of the lens.

PHYSICAL CHARACTERISTICS:

Lens / Haptic Material	Modelo SIP: Foldable hydrophilic acrylic copolymer containing a UV absorber. Modelo SIPY: Foldable hydrophilic acrylic copolymer containing a UV absorber and a blue light filtering chromophore.
Material Characteristics	Index of Refraction when wet at 20° C: 1.460 <p>Index of Refraction when in the eye at 35° C: 1.459</p>
Spectral Transmittance	Modelo SIP: 10% transmittance at 365 nanometers for the “20 diopter IOL.” <p>Modelo SIPY: 10% transmittance at 374 nanometers for +20.0 diopter IOL. See FIG. 8 with chart’s X-value = Wavelength (nm) and Y-value = % Transmittance</p>
Optic Type	Aspheric aberration free
Powers	-10.0 to +40.0 Diopters
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 10.5 mm to 11 mm across the dioptr range; Haptic Angle: 0°
Lens Orientation	The anterior side of the lens is facing up. There are various surgical procedures that can be utilized. As illustrated in FIG. 7, when the orientation features are top right (A) and bottom left (B). (FIG. 7).

When the “click-lock” mechanism is engaged, the lens is securely in place. The system is ready for injection. For detailed description, refer to the instructions for use supplied with the compatible injector.

11.Check that the lens is positioned centrally in the loading chamber and close the wings of the loading chamber (FIG. 6). When the “click-lock” mechanism is engaged, the lens is securely in place. The system is ready for injection.

12.Press the injector plunger forward to push the lens. Implant the lens with the bevel of the injector tip facing down to the patient’s eye. The anterior side of the lens is facing up towards the cornea while the viscoelastic features are top right (A) and bottom left (B) (FIG. 7).

Only viscoelastic solution of low to medium viscosity should be used. Only the injection system provided with the preloaded Synthesis lens should be used.

There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient.

SYMBOLS USED ON LABELING:

Symbol	Description
D	Diopter
ØB	Body Diameter
ØT	Overall Diameter

INDICATIONS:
The Synthesis posterior chamber lens is indicated for primary implantations for the visual correction of aphakia in capsular lenses. The lens is designed for implantation in the capsular bag.

CONTRA-INDICATIONS AND WARNINGS:
Physicians considering lens implantation in any of the following circumstances (but not limited to) must weigh the potential risk/benefit ratio: recurrent severe anterior or posterior segment inflammation or uveitis; patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases; surgical difficulties increasing the potential for complications (such as but not limited to: persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss); distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL would result in damage to the lens; patients who have the endothelium during implantation; patients in whom neither the posterior capsule nor zonule are intact enough to provide support; suspected microbial infection.
As with any surgical procedure, there is risk involved. Potential complications of cataract surgery may include, but are not limited to the following: inflammation (iridocyclitis), pupillary membrane, vitreous inflammation, CME, TASS, ...), corneal endothelial damage, endothelialmitis, retinal detachment, corneal edema, pupillary block, iris prolapse, hyppopyon, glaucoma, wound leak, hyphema, iris atrophy, vitreous prolapse, IOL decentration, or tilt. IOL silicone adhesives, IOL optic deposits or opacities, and secondary surgical intervention (such as but not limited to wound leak repair, vitreous aspiration, iridectomy, IOL repositioning or replacement and retinal detachment repair).
The IOL must be implanted in accordance to the following instructions for use. Improper use may impose risk to the patient’s health.
Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage. Do not attempt to modify or alter the device or any of its components.
Do not use lens after expiration date.
Do not reuse the lens and the preloaded system. Reuse may impose serious risk to the patient’s health.
Do not resterilize as this can produce undesirable side effects.

- PRECAUTIONS:**
- Store at room temperature. Avoid temperatures below 5° C and above 45° C. Do not freeze.
 - The IOL model, its power and expiration date should be verified before opening the protective packaging and before opening the individual sterile pouch.
 - Use correct injector compatible.
 - Sterility of the IOL and of the injector parts is guaranteed only if the individual sterile pouches have not been opened or damaged.
 - IOLs must be carefully rinsed prior to implantation. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution.
 - Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens after implantation in the capsular bag.
 - If a YAG laser posterior capsulotomy is performed, assure that the laser beam is focused slightly behind the posterior capsule.

PACKAGING / STERILIZATION:
The Preloaded Synthesis IOL is individually packaged in a pouch and a vial containing balanced salt solution that should be used under sterile conditions. The Synthesis IOL is placed in a loading chamber, fixed by a lens holder, all immersed in BSS. A patient card and self-adhesive labels are supplied in the carton to provide traceability of the IOL. The Preloaded Synthesis IOL is steam sterilized.

DIRECTIONS FOR USE:

1. The physician should determine preoperatively the power of the lens to be implanted, using calculation methods. Recommended A-constants listed on the lens carton are guidelines only.

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