Aspheric preloaded hydrophobic lens (EN) – [Instruction For Use \(en\) - LuxGood™](#)

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CE0197

INSTRUCTION FOR USE	
Date of the last revision of IFU: 2023-03-07	Aspec-e19462-04

1. Product to which these instructions for use apply


These instructions for use apply to the following products:

Brand Name	Models	Performance characteristics	Basic UDI-DI
LUXGOOD	GOOD (LUXGOOD™ Crystal)	Aspheric preloaded hydrophobic lens	3665493LUXGS8K
	YGOOD (LUXGOOD™)		

2. Description

The LuxGood™ intraocular lens (IOL) is a single-piece ultraviolet (UV) absorbing posterior chamber IOL developed to replace the natural crystalline lens in aphakic adult patients. The LuxGood™ IOL has aspheric optics. The lens is provided sterile and is designed to be folded prior to insertion in the eye and implantation in the capsular bag. The preloaded intraocular LuxGood™ lens is supplied with a delivery system to provide a touch-free method of delivery of the lens into the capsular bag. The LuxGood™ IOL is supplied in 2 different models – Models GOOD or YGOOD. See label on the cardboard box for type of the lens.

3. Physical Characteristics

Lens / Haptic Material	Model GOOD: Foldable hydrophobic acrylic copolymer containing a UV absorber. Model YGOOD: Foldable hydrophobic acrylic copolymer containing a UV absorber and a blue light filtering chromophore.
Material Characteristics	LuxGood™ hydrophobic acrylic IOLs are made from highly purified hydrophobic acrylic copolymer composed of 56% ethylene glycol phenyl ether acrylate (EGPEA), 21% butyl acrylate (BA), 16.6% methylstyrene (MeS) and 2% hydroxyethyl methacrylate (HEMA), cross-linked with ethylene glycol dimethacrylate (EGDMA) including a benzotriazole UV absorber. Yellow dye is added to the material of Model YGOOD. Index of Refraction when wet at 21°C: 1.544 Index of Refraction when in the eye at 35°C: 1.540
Spectral Transmittance	Model GOOD: 10% transmittance at 393.5 nanometers for +20.00 diopter IOL. Model YGOOD: 10% transmittance at 396.7 nanometers for +20.00 diopter IOL.  with chart's X-value = Wavelength (nm) and Y-value = % Transmittance
Optic Type	Models LuxGood™ Crystal and LuxGood™: Aspheric aberration free
Powers	+00.00 to +34.00 Diopters
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 11 mm; Haptic Angle: 0°
Lens Orientation	The IOL is to be implanted with the anterior side of the lens facing up towards the anterior side of the eye. The orientation of the IOL can be verified by visual inspection of the lens. As illustrated in FIG.5, when the orientation features are top right (A) and bottom left (B), you are facing the anterior side of the lens.

4. Indications

Intended purpose:

LuxGood™ IOLs are intended to replace the human crystalline lens.

Indication:

The LuxGood™ posterior chamber lens is indicated for primary implantation for the visual correction of aphakia in patients in whom the cataractous lens has been removed. The lens is designed for implantation in the capsular bag.

Patient target group:

Aphakic adult patients.

Intended users:

IOLs must be handled by health professionals and implanted by physicians.

Mode of action:

LuxGood™ IOL positioned in the capsular bag acts as a refractive medium to correct aphakia.

Lifetime:

The LuxGood™ IOL estimated lifetime is 20 years.

5. Clinical benefits

The clinical benefit of the implantation of an IOL for cataract patients is the prevention of blindness.

6. Safety and clinical performance

The summary of safety and clinical performance (SSCP) is published in the European Database on Medical Devices (Eudamed).

The URL of the public website of Eudamed is: <https://ec.europa.eu/tools/eudamed>.

The SSCP is linked to the Basic UDI-DI listed in section 1.

7. Packaging/Sterilization

The LuxGood™ IOL is individually package preloaded in a delivery system in a blister. All models are supplied dry and should be opened under sterile conditions. A patient card and self-adhesive labels are supplied in the carton to provide traceability of the IOL. All models are sterilized by Ethylene Oxide.

8. Precautions for use and storage

- Store at room temperature. Avoid temperatures below 5°C and above 45°C. Do not freeze.
- The LuxGood™ IOL, the balanced salt solution and the viscoelastic solution used should be stored for at least 30 min at minimum 20°C before surgery.
- A high level of surgical skill is required to implant intraocular lenses. Surgeon should be fully trained before attempting to implant intraocular lenses.
- The IOL model, its power and expiration date should be verified before opening the protective packaging and before opening the individual sterile blister.
- Sterility of the device is guaranteed only if the individual sterile blister has not been opened or damaged.
- Handle the lens carefully to avoid surface or haptics damages.
- The IOL must be carefully rinsed prior to implantation. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution.
- Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens after implantation in the capsular bag.
- Care should be taken to verify that the lens is fully implanted into the capsular bag.
- If a YAG laser posterior capsulotomy is performed, assure that the laser beam is focused slightly behind the posterior capsule.

Disposal

Disposal of the injection system and discarded IOLs (used or unused) are classified as medical (clinical) waste that harbors a potential infection or microbial hazard and must be disposed of accordingly.

Disposal should be done in compliance with the local regulation of waste from health-care activities characterized by risk of infection.

9. Directions for use

1. The physician should determine preoperatively the power of the lens to be implanted, using calculation methods. The recommended A-constants listed on the lens carton are guidelines only.
2. Prior to implanting, examine the lens package for model, power, expiration date and proper configuration.
3. Open the carton and inspect the sterilized blister containing the lens and the injection system. Make sure it is not damaged. Verify that the information on the blister label is consistent with the information on the outer box and the self-adhesive labels. Peel open the blister and transfer the device onto the sterile field.
4. Inspect the device and make sure it is intact. The lens should be rinsed thoroughly by flushing balanced salt solution through the tip of the cartridge (FIG. 1) and from the end of the loading chamber (FIG. 2).

Please refer to the enclosed figures to fully prepare the device before implantation, taking care not to touch the IOL:

FIG. 3: Apply viscoelastic only in the tip of the cartridge. Wait for 30 seconds.


FIG. 4: Close the wings of the loading chamber. When the "click-lock" mechanism is engaged, the lens is securely in place. The lens is ready for implantation.

5. Inject the lens in the shortest possible time after fully preparing the lens.

6. Insert the cartridge tip in the incision with the bevel of the injector tip facing down to the patient's eye. Proceed by pushing the plunger forward smoothly but continuously until the lens is delivered into the capsular bag.

7. Discard the injection device (see §8 Disposal).

8. The lens is to be implanted with the anterior side of the lens facing up towards the anterior side of the eye. The orientation of the IOL can be verified by visual inspection of the haptics. As illustrated in FIG. 5, when the haptic insertions are top right (A) and bottom left (B), you are looking at the anterior side of the lens.

There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient.

10. Contra-indications and Warnings

Physicians considering lens implantation in any of the following circumstances (but not limited to) must weigh the potential risk/benefit ratio: recurrent severe anterior or posterior segment inflammation or uveitis; patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases; surgical difficulties increasing the potential for complications (such as but not limited to: persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss); distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible; circumstances that would result in damage to the endothelium during implantation; patients in whom neither the posterior capsule nor zonule are intact enough to provide support; suspected microbial infection.

Warnings

The IOL must be implanted in accordance with the instructions for use. Improper use may impose risk to the patient's health.

Do not use if product sterility or quality is thought to be compromised due to damaged packaging or if product was unintentionally opened before use.

Do not attempt to modify or alter the device or any of its components.

Do not use after expiration date.

Do not use without viscoelastic solution.

Do not reuse. Reuse may impose serious risk to the patient's health.

Do not resterilize as this can produce undesirable side effects.

Do not implant the lens in the anterior chamber or in the ciliary sulcus.

Only viscoelastic solution of low to medium viscosity should be used.

Only the injection system provided with the preloaded LuxGood™ lens should be used.

11. Possible undesirable side effect

As with any surgical procedure, there is risk involved. Potential complications of cataract surgery may include, but are not limited to the following: secondary cataract, inflammation (iridocyclitis, pupillary/cyclitic membrane, vitreous inflammation, CME, TASS ...), elevated IOP, anterior capsule opacification (fibrosis, phimosis), corneal endothelial damage, endophthalmitis, retinal detachment, corneal edema, pupillary block, iris prolapse, hypopyon, glaucoma, wound leak, hyphema, iris atrophy, vitreous prolapse, IOL decentration or tilt, IOL optic deposits or opacities, glistenings, residual refractive error, visual symptoms and secondary surgical intervention (such as but not limited to wound leak repair, vitreous aspiration, iridectomy, IOL repositioning or replacement and retinal detachment repair).

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

12. Implant card and patient information

The implant card included in the package is to be completed and given to the patient together with instructions to keep this card as a permanent record of the implant and to show the card to any eye care professional consulted in the future.

Patient information is made available on the internet. The link to access the information is printed on the implant card.

13. Information to be conveyed to the patient

- Patients should be advised that unexpected outcomes could lead to continued spectacle dependence.

14. Symbol used on labelling

	Manufacturer		Do not re-use
	Date of manufacture		Consult instructions for use or consult electronic instructions for use
	Use-by date		Caution
	Batch code		Patient name
	Catalogue number		Health care center or doctor
	Serial number		Date
	Sterilized using ethylene oxide		Medical device
	Do not re-sterilize		Unique device identifier
	Do not use if package is damaged and consult instructions for use	\emptyset_B	Diameter of the body
	Single sterile barrier system with protective packaging outside	\emptyset_T	Overall diameter
	Lower limit of temperature	D	Diopter
	Right eye		Left eye

 Cutting Edge S.A.S.
770 rue Alfred Nobel
Immeuble Le Nobel
34000 Montpellier
FRANCE

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