Trifocal Toric preloaded hydrophobic lens (EN) – [Instruction for Use \(en\) – LuxLife Toric](#)
Lentille hydrophobe trifocale torique préchargée (FR) – [Notice d'utilisation \(fr\) – LuxLife Toric](#)

 Cutting Edge S.A.S.
770 rue Alfred Nobel
Immeuble Le Nobel
34000 Montpellier
FRANCE

CE 0197

INSTRUCTION FOR USE

Date of the last revision of IFU: 2023-07-20

Ref: Aspec-e19851-00

1. Product to which these instructions for use apply


These instructions for use apply to the following products:

Brand Name	Models	Performance characteristics	Basic UDI-DI
LUXLIFE Toric	LIFET(LuxLife Toric Crystal)	Trifocal Toric preloaded hydrophobic lens	3665493LUXLIFEFD
	YLIFET (LuxLife Toric yellow)		

2. Description

The LuxLife Toric intraocular lens (IOL) is a single-piece ultraviolet (UV) absorbing posterior chamber IOL developed to replace the natural crystalline lens in aphakic patients with pre-existing corneal astigmatism. The LuxLife Toric IOL has trifocal biconvex aspheric and toric optics with refractive optic structure providing increased depth of focus. The near add power of +3.4 diopters and intermediate add power of +2.2 diopters provide increased spectacle independence for intermediate and near distances. The LuxLife Toric™ lens incorporates a pure refractive optics (PRO) technology, across its entire optical diameter with no diffractive rings. The lens is provided sterile and is designed to be folded prior to insertion in the eye and implantation in the capsular bag. The preloaded intraocular LuxLife Toric lens is supplied with an Accuject™ Dual Delivery System to provide a touch-free method of delivery of the lens into the capsular bag. The LuxLife Toric IOL is supplied in 2 different models – Models LIFET or YLIFET. See label on the cardboard box for type of the lens.

3. Physical Characteristics

Lens / Haptic Material	Model LIFET: Foldable hydrophobic acrylic copolymer containing a UV absorber. Model YLIFET: Foldable hydrophobic acrylic copolymer containing a UV absorber and a blue light filtering chromophore.
Material Characteristics	LuxLife Toric hydrophobic acrylic IOLs are made from highly purified hydrophobic acrylic copolymer composed of 56% ethylene glycol phenyl ether acrylate (EGPEA), 21% butyl acrylate (BA), 16.6% methylstyrene (MeS) and 2% hydroxyethyl methacrylate (HEMA), cross-linked with ethylene glycol dimethacrylate (EGDMA) including a benzotriazole UV absorber. Yellow dye is added to the material of Model YLIFET. Index of Refraction when wet at 21°C: 1.544 Index of Refraction when in the eye at 35°C: 1.540
Spectral Transmittance	Model LIFET: 10% transmittance at 393.5 nanometers for +20.0 diopter IOL. Model YLIFET: 10% transmittance at 397 nanometers for +20.0 diopter IOL.  <p style="text-align: center;">X-value = Wavelength (nm) and Y-value = % Transmittance</p>
Optic Type	Refractive aspheric trifocal toric Intermediate Add power: +2.2 diopters Near Add power: +3.4 diopters

Powers (SE)	+06.00 to +34.00 Diopters
Cylindrical powers	+0.75 to +6.00 Diopters
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 11 mm; Haptic Angle: 0°
Lens Orientation	The IOL is to be implanted with the anterior side of the lens facing up towards the anterior side of the eye. The orientation of the IOL can be verified by visual inspection of the lens. As illustrated in FIG.4, when the orientation features are top right (A) and bottom left (B), you are facing the anterior side of the lens.
MTF curve	<p>Through focus MTF curve of the LuxLife Toric IOL at 50-line pairs / mm using the ISO model eye and an aperture of 2/3/5 mm. This graph determines the optical depth of focus of the lens with higher values representing the higher performances for wavelength ranging from 300 nm to 1100 nm.</p> <p style="text-align: center;">MTF Through Focus Response @ 50lp/mm for 2,0mm, 3,0mm and 5mm aperture</p> <p style="text-align: center;">Legend: 2.0 mm (orange), 3.0 mm (blue), 5.0 mm (grey)</p>

4. Indications

Intended purpose:

This LuxLife Toric IOL is intended for implantation in the capsular bag to replace the human crystalline lens.

Indication:

The LuxLife Toric posterior chamber lens is indicated for primary implantation for the visual correction of aphakia and pre-existing corneal astigmatism in patients in whom the cataractous lens has been removed and who desire near, intermediate and distance vision with increased spectacle independence.

Patient target group:

Aphakic adult patients with pre-existing regular corneal astigmatism.

Intended users:

IOLs must be handled by health professionals and implanted by physicians.

Mode of action:

LuxLife Toric IOL positioned in the capsular bag acts as a refractive medium to correct aphakia and pre-existing corneal astigmatism. This IOL has a proprietary biconvex refractive trifocal optic that provides increased depth of focus.

LuxLife Toric IOL has a posterior toric optical surface with axis marks that identify the flat meridian of the cylinder (plus cylinder axis) (FIG.4). Alignment of the axis marks with postoperative steep corneal meridian (intended axis of placement) ensures correct placement and optimal results.

5. Clinical benefits

The clinical benefit of the implantation of an IOL for cataract patients is the prevention of blindness.

The LuxLife Toric IOL provides functional near, intermediate, and far vision, and decreases patients' dependence on glasses.

6. Safety and clinical performance

The summary of safety and clinical performance (SSCP) is published in the European Database on Medical Devices (Eudamed).

The URL of the public website of Eudamed is: <https://ec.europa.eu/tools/eudamed>.

The SSCP is linked to the Basic UDI-DI listed in section 1.

7. Summary of clinical data

The premarket clinical study “Multicentric Clinical Study to Determine Safety and Efficacy of a Hydrophobic Acrylic Trifocal Intraocular Lens (IOL)” (NCT04465344) was performed in the EU at up to 7 investigational sites. It is a multicentric, non-randomized, controlled, open-label ethics-approved trial with bilateral implantation of 188 subjects in total (376 eyes - 125 subjects implanted with the investigational LuxLife IOL and 63 subjects implanted with the CE approved LuxGood monofocal IOL) with a follow-up of 6 months after cataract surgery.

The objective was to demonstrate safety and performance of the device, comparing SPE rates of the investigational device to the results of historical data of EN ISO 11979-7:2018 and demonstrating non inferiority of LuxLife compared to the monofocal control lens in terms of best corrected distance visual acuity.

Results demonstrated safety of the LuxLife lens. Clinical performance was comparable between LuxLife and LuxGood for distance vision attesting for the lens ability to compensate for aphakia. A statistically significant increase between pre- and postoperative examination on monocular Distance Corrected Near Visual Acuity and Distance Corrected Intermediate Visual Acuity demonstrated the ability of the LuxLife IOL to provide adequate intermediate and near visual acuity with low spectacle dependency (shown by patient's reported outcomes). Contrast sensitivity and photic phenomena outcomes were worse with LuxLife but in a usual and expected range.

For further information on available clinical data, please refer to the SSCP available online (chapter 6).

8. Packaging/Sterilization

The LuxLife Toric IOL is individually package preloaded in an Accuject™ Dual delivery system in a blister. All models are supplied dry and should be opened under sterile conditions.

A patient card and self-adhesive labels are supplied in the carton to provide traceability of the IOL.

The LuxLife Toric IOL is labeled with spherical equivalent and cylindrical powers. Spherical equivalent corresponds to the power obtained with most calculation formulas.

All models are sterilized by Ethylene Oxide.

9. Precautions for use and storage

- Store at room temperature. Avoid temperatures below 5°C and above 45°C. Do not freeze.
- The LuxLife Toric IOL, the balanced salt solution and the viscoelastic solution used should be stored for at least 30 min at minimum 20°C before surgery.
- A high level of surgical skill is required to implant intraocular lenses. Surgeon should be fully trained before attempting to implant intraocular lenses.
- The IOL model, its powers and expiration date should be verified before opening the protective packaging and before opening the individual sterile blister.
- Sterility of the device is guaranteed only if the individual sterile blister has not been opened or damaged.
- Handle the lens carefully to avoid surface or haptics damages.
- The IOL must be carefully rinsed prior to implantation. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution.
- Care should be taken to verify that the lens is fully implanted into the capsular bag.
- If a YAG laser posterior capsulotomy is performed, assure that the laser beam is focused slightly behind the posterior capsule.
- Customization of surgically induced astigmatism (SIA) is recommended based upon individual surgical technique and past results since size and location of the surgical incision may affect amount and axis of postoperative corneal astigmatism.
- Consideration of posterior corneal astigmatism is recommended for optimal results.
- Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens after implantation in the capsular bag, not to overinflate the capsular bag and to ensure that the wound is watertight at the end of the case. Residual viscoelastic / overinflation / wound leak may allow the lens to rotate away from the intended axis of placement.
- Special care should be taken to ensure proper orientation of the toric lens at the intended axis of placement at the end of the surgical case.

Disposal

Disposal of the injection system and discarded IOLs (used or unused) are classified as medical (clinical) waste that harbors a potential infection or microbial hazard and must be disposed of accordingly.


Disposal should be done in compliance with the local regulation of waste from health-care activities characterized by risk of infection.

10. Directions for use

1. Accurate biometry is essential to successful visual outcomes. The physician should determine preoperatively the spherical equivalent and cylindrical powers of the lens to be implanted, using

calculation methods and accurate biometry measurements depending on the surgeon's experience, preferences and/or with the assistance of a Toric Calculator. The recommended A-constants listed on the lens carton are guidelines only.

2. Prior to surgery, using the physician's preferred method, the operated eye should be referenced to provide a clear indication of the intended axis of placement.
3. Prior to implanting, examine the lens package for model, powers (spherical equivalent and cylinder), expiration date and proper configuration.
4. Open the carton and inspect the sterilized blister containing the lens and the injection system. Make sure it is not damaged. Verify that the information on the blister label is consistent with the information on the outer box and the self-adhesive labels. Peel open the blister and transfer the device onto the sterile field.
5. Inspect the device and make sure it is intact. The lens should be rinsed thoroughly by flushing balanced salt solution through the tip of the cartridge (FIG. 1a) and from the end of the loading chamber (FIG. 1b). Please refer to the enclosed figures to fully prepare the device before implantation, taking care not to touch the IOL:
Apply viscoelastic in the tip of the cartridge (FIG. 1a). Wait for 30 seconds.
Close the wings of the loading chamber (FIG. 2). When the "click-lock" mechanism is engaged, the lens is securely in place. The lens is ready for implantation.
6. Press the injector plunger forward until the plunger stops in the middle of the injector housing (FIG. 3). The lens should be positioned in the housing of the injector until implantation and should be injected in the shortest possible time after closing the loading chamber.
7. Insert the cartridge tip in the incision with the bevel of the injector tip facing down to the patient's eye.
In push mode: open the wings and inject by pushing the plunger forward smoothly but continuously.
In screw mode: inject by rotating the screw mechanism clockwise to advance the lens forward until the lens is delivered into the capsular bag.
8. Discard the injection device.
9. The lens is to be implanted with the anterior side of the lens facing up towards the anterior side of the eye. The orientation of the IOL can be verified by visual inspection of the haptics. As illustrated in FIG. 4, when the haptic insertions are top right (A) and bottom left (B), you are looking at the anterior side of the lens.
10. Align LuxLife Toric IOL's axis marks with the referenced intended axis of placement.

There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient.

11. Contra-indications and Warnings

Physicians considering lens implantation in any of the following circumstances (but not limited to) must weigh the potential risk/benefit ratio: recurrent severe anterior or posterior segment inflammation or uveitis; significant corneal irregularities or dystrophy; amblyopia; low potential for good vision; iris or pupil abnormalities, glaucoma, patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases; surgical difficulties increasing the potential for complications (such as but not limited to: persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss); distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible; circumstances that would result in damage to the endothelium during implantation; patients in whom neither the posterior capsule nor zonule are intact enough to provide support and stability; suspected microbial infection.

Toric IOLs are not recommended with pre-existing irregular astigmatism, radial tears on the capsulotomy, zonular damage or when the need for a large capsulotomy is anticipated.

Patients with preoperative problems such as corneal endothelial disease, abnormal cornea, macular degeneration, retinal degeneration, glaucoma, amblyopia, and chronic drug miosis may not achieve the visual acuity of patients without such problems. Identification, prevention, and treatment of ocular surface disease can help improve patient satisfaction with multifocal IOLs.

The surgeon must target emmetropia to achieve optimal visual performance.

Patients with predicted postoperative astigmatism > 1.0 D may not achieve optimal visual outcomes.

Care should be taken to achieve IOL centration, as lens decentration may result in a patient experiencing visual disturbances under certain lighting conditions.

Some visual effects may be expected due to the superposition of focused and unfocused multiple images.

These may include some perceptions of halos or radial lines around point sources of light under nighttime conditions. A reduction in contrast sensitivity as compared to a monofocal IOL may be experienced by some patients. On rare occasions, these visual effects may be significant enough that the patient may request removal of the IOL.

When performing refraction, interpret results with caution when using autorefractors or wavefront aberrometers. Confirmation of refraction with maximum plus manifest refraction technique is recommended.

Warnings

The IOL must be implanted in accordance with the following instructions for use. Improper use may impose risk to the patient's health.

Do not use if product sterility or quality is thought to be compromised due to damaged packaging or if product was unintentionally opened before use.

Do not attempt to modify or alter the device or any of its components.

Do not use after expiration date.

Do not use without viscoelastic solution.

Do not reuse. Reuse may impose serious risk to the patient's health.

Do not resterilize as this can produce undesirable side effects.

Do not implant the lens in the anterior chamber or in the ciliary sulcus.

Rotation of LuxLife Toric lens away from its intended axis can reduce its astigmatic correction.

Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, repositioning should occur as early as possible prior to lens encapsulation.

Do not use refractive data for calculation since the presence of lenticular astigmatism in the natural lens to be removed may influence results.

Only viscoelastic solution of low to medium viscosity should be used.

Only the injection system provided with the preloaded LuxLife Toric lens should be used.

12. Possible undesirable side effect

As with any surgical procedure, there is risk involved. Potential complications of cataract surgery may include, but are not limited to the following: secondary cataract, inflammation (iridocyclitis, pupillary/cyclitic membrane, vitreous inflammation, CME, TASS...), anterior capsule opacification (fibrosis, phimosis), corneal endothelial damage, endophthalmitis, retinal detachment, corneal edema, pupillary block, iris prolapse, hypopyon, glaucoma, wound leak, hyphema, iris atrophy, vitreous prolapse, IOL decentration or tilt, dysphotopsia, IOL optic deposits or opacities, glistenings, residual refractive error, visual symptoms and secondary surgical intervention (such as but not limited to wound leak repair, vitreous aspiration, iridectomy, IOL repositioning or replacement and retinal detachment repair).

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

13. Implant card and patient information

The implant card included in the package is to be completed and given to the patient together with instructions to keep this card as a permanent record of the implant and to show the card to any eye care professional consulted in the future.


Patient information is made available on the internet. The link to access the information is printed on the implant card.

14. Information to be conveyed to the patient

- Some patients may experience reduced contrast sensitivity as compared to monofocal IOLs.
- Some patients may experience visual effects with the Trifocal IOLs because of the super-positioning of focused and unfocused images. Visual effects may include the perception of halos or radial lines around point light sources under low illumination conditions.
- Patients should be advised that unexpected outcomes could lead to continued spectacle dependence.

15. Symbol used on labelling

	Manufacturer		Do not re-use
	Date of manufacture		Consult instructions for use or consult electronic instructions for use

	Use-by date		Caution
	Batch code		Patient name
	Catalogue number		Health care center or doctor
	Serial number		Date
	Sterilized using ethylene oxide		Medical device
	Do not resterilize		Unique device identifier
	Do not use if package is damaged and consult instructions for use	\varnothing_B	Diameter of the body
	Single sterile barrier system with protective packaging outside	\varnothing_T	Overall diameter
	Lower limit of temperature	D	Diopter
	Right eye		Left eye

 Cutting Edge S.A.S.
770 rue Alfred Nobel
Immeuble Le Nobel
34000 Montpellier
FRANCE

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