




CUTTING EDGE
sense of innovation



Synthesis[®] Toric




High
precision



CUTTING EDGE
sense of innovation

INNOVATION IN ACTION



2017

Synthesis®
Acrylic Monofocal


Genesis®
Hydrophobic Monofocal



2018

Synthesis®
EDOF

Synthesis[®] Toric



2019

The **Synthesis[®] Toric** is engineered by a Research and Development team experienced and dedicated to the design and manufacture of high-quality intraocular lenses.


The **Synthesis[®] Monofocal**, **Synthesis[®] \oplus EDOF** and **Synthesis[®] Toric** comprise of a unique family of micro-incision IOLs designed for the health care professional and to make a contribution to patient well-being.

Managing astigmatism is made easier by combining the strengths of: a proven rotational stability IOL platform, a large combination of sphere/cylinder options and the latest intuitive **Synthesis[®] Toric** online calculator.

ASTIGMATISM PATIENT CARE

- **Astigmatism needs to be treated**
to reduce the dependency on glasses for distance vision.⁽¹⁾
- **78% of the cataract patient population present with corneal astigmatism $\geq 0,5$ D.**⁽²⁾

Distribution of preoperative astigmatism
in patients undergoing cataract surgery (N=110,468)⁽²⁾



TORIC IOLS CAN OFFER AN EFFECTIVE SOLUTION.⁽⁴⁾

THE ADOPTION AND USAGE OF TORIC IOLS REMAIN RELATIVELY LOW DEPENDING ON THE COUNTRY, RANGING FROM (4 to 9%).⁽³⁾

Takes into consideration
all sources of astigmatism ⁽⁵⁾

Utilises the latest generation
online calculator ⁽⁶⁾

Offers a large
combination of
sphere/cylinder

Ensures
rotational
stability ⁽⁷⁾

Minimises
surgically
induced
astigmatism ⁽⁸⁾

High
precision

NEXT GENERATION ONLINE CALCULATOR

- Takes into account posterior corneal astigmatism.
 - Ignoring posterior corneal astigmatism may lead to the incorrect estimation of the total corneal astigmatism.⁽⁵⁾

CUTTING EDGE CUTTING-EDGE Toric Calculator English Log in

Information

Surgeon's name: James Patient's ID: 123456 Date of birth: 10.10.1952 67 y.o. Gender: ☒ M ☐ F

IOL Model: Synthesis Toric Eye: ☒ OD ☐ OS Include posterior astigmatism: ☒ Yes ☐ No

Flat K: 42.51 D @ 123° K index: 1.3375 Steep K: 44.12 D @ 33° SIA: 0.1 D Incision: 12°

AL: 23.41 mm ACD: 3.1 mm A-constant: 118.85 Refractive target: -0.5 D

Ensure that the data entered are correct. Cutting-Edge does not send or save information for patient identification. By clicking on the 'Accept and calculate' button, you agree to the **TERMS AND CONDITION OF USE**.

Confirm and Calculate Clear all

Expected Results

Cylinder 1.0 D Cylinder 1.5 D Cylinder 2.25 D

IOL power (SE) (D)	SPK/ST	IOL orientation	Residual astigmatism
+22.25	-0.46	34°	0.52D @ 34°
+22.75	-0.82	34°	0.52D @ 34°
+23.25	-1.17	34°	0.52D @ 34°

Pulesance IOL recommandé (SE): 22.75 0.52D @ 34°

Print Quotation request

CUTTING EDGE Instruction for use General conditions of use


toriccalc.cutting-edge.fr



CLINICAL RESULTS


ON ROTATIONAL STABILITY

- **Primary objective:**
Synthesis 11 mm platform rotational stability analysis with 4 months follow-up.⁽⁹⁾
- **Prospective monocentric clinical study.**
- **Clinically proven rotational stability:**
93% of patients with less than 5° rotation from surgery day (D0)
over a 4 month follow-up period.
Average rotation of 2.93° at 4 months.



2,93°⁺ or - 1,99°


AVERAGE
ROTATION
BETWEEN
D0 AND
4 MONTHS




LARGE RANGE

- **Sphere**
from +5,00D to +32,00D
- **Cylinder**
from +0.75D, +1.0D, +1.50D,
+2.25D, +3.00D, +3.75D, +4.50D,
+5.25D, 6.00D
- **MiCS** 1,6 mm
- **Synthesis® Toric** can treat corneal astigmatism from 0.50D.

PRELOADED INJECTION SYSTEM



Synthesis® Toric

		Brands	References	Wound-assisted recommended incision	In-the wound recommended incision	In-the-bag recommended incision
Synthesis® Toric	Preloaded	Medicel Accuject™ PRO 1.6	LP604550P (box of 10)	1.6 mm	1.8 mm	2.0 mm

TECHNICAL SHEET

Synthesis[®] Toric



References	PRELOADED
Clear implant	SIPT

DESIGNATION	SPECIFICATIONS
Optic diameter	6 mm across the power range
Total diameter	11 mm across the power range
Design	4 point-fixation, C-Edge+ ≤ 5 microns
Optical	Aspheric and toric without aberration
Angulation	0°, posterior-shift feature
Dioptric power (Sphere)	1.0D increment between +5.00D to +10.00D 0.5D increment between +10.50D to +30.00D 1.0D increment between +31.00D to +32.00D
Dioptric power (Cylinder)	0,75D/1,0D/1,5D/2,25D/3,0D/3,75D/4,5D/5,25D/6,0D
Biometry available on the IOLcon website https://iolcon.org/	Optical biometry: SRK/T: A=118,85 Haigis: a0=1,223 (with a1=0,4 and a2=0,1) Hoffer-Q: pACD=5,46 Holladay 1: SF=1,706
Incision size	From 1.6 mm

toriccalc.cutting-edge.fr

This product is guaranteed latex free and phthalates free.

Synthesis[®] Toric

HIGH PRECISION

Next generation
online calculator




MICS
for surgically induced astigmatism



Large combination
of Sphere/Cylinder



Clinically proven
rotational stability



CUTTING EDGE

LOADING GUIDE

Synthesis[®] Toric

PRELOADED


References

PRELOADED

Synthesis Toric Clear

SIPT

- ① Open the container; remove the loading cartridge that contains the IOL from the container by pulling on the armature taking care not to remove the protective clip.
- ② Insert the loading cartridge into the injector body until it locks into position.
- ③ Carefully rinse the lens with balanced salt solution.




- ④ Through the filling hole, fill the cartridge tunnel with viscoelastic.
- ⑤ Fill a small amount of visco behind the lens at the rear of the loading chamber.
- ⑥ Push the plunger until the silicone tip approaches the rear of the loading chamber.




- ⑦ Carefully remove the clip from the loading chamber:
 - Grab the armature with two fingers while rotating the bottom hinge upwards with the thumb
 - Lift vertically.



- ⑧ Check that the lens is correctly positioned and close the wings of the loading chamber. When loaded correctly, you will hear a "click" sound.
- ⑨ Push the plunger to position the implant in the cartridge tunnel. The system is ready for injection.



- ⑩ Once implanted, the orientation features of the Synthesis[®] Toric's implant must be oriented at the top right (A) and at the bottom left (B).



An innovative family of IOLs



PRELOADED | YELLOW | CLEAR | MONOFOCAL | EDOF | TORIC

CUTTING EDGE

770 rue Alfred Nobel - Immeuble le Nobel - 34000 Montpellier - France

Tél. : +33 (0)5 62 24 65 03 - Fax +33 (0)5 62 24 65 05

KESTREL OPHTHALMICS LTD

Unit 9 Cabot Business Village, Holyrood Close, Poole, Dorset, BH17 7BA - United Kingdom

Tel. +44 (0) 1202 696963

CUTTING EDGE SURGICAL

70 izq Local Avinguda Manuel Girona - 08860 Casteldefels - España

Tel. +34 931 293 849

www.cutting-edge.eu - info@cutting-edge.eu



Bibliography >

- (1) Wilkins M, Allan B, Rubin G. Spectacle use after routine cataract surgery. Br J Ophthalmol 2009;93:1307-1312
- (2) Day AC, Dhariwal M, Keith MS, et al. Distribution of preoperative and postoperative astigmatism in a large population of patients undergoing cataract surgery in the UK. British Journal of Ophthalmology. Published Online First: 06 September 2018. doi: 10.1136/bjophthalmol-2018-312025
- (3) Freeman B. Market scope. 2016 IOL Report: A Global Market Analysis for 2015 to 2021.
- (4) Kessel L, Andresen J, Tendal B et al. Toric Intraocular Lenses in the Correction of Astigmatism During Cataract Surgery. A Systematic Review and Meta-analysis. Ophthalmology. 2016. 123(2):275-86
- (5) Koch D, Ali S, Weikert M et al. Contribution of posterior corneal astigmatism to total corneal astigmatism. J Cataract Refract Surg 2012; 38:2080-2087
- (6) Levy P. «Pré requis , implants toriques calcul du cylindre» . Amzallag T, Rozot P. Chirurgie de la cataracte. 2018, Elsevier Masson SAS. Pages 19 à 24 .
- (7) Kaur M, Shaikh F, Falera R. et al. Optimizing outcomes with toric intraocular lenses. Indian J Ophthalmol 2017;65:1301-13
- (8) Tetikoglu M, Yeter C, Helvacioğlu F et al. Effect of Corneal Incision Enlargement on Surgically Induced Astigmatism in Bifocal Microincision Cataract Surgery. Turk J Ophthalmol 2016. 46; 3: 99-103
- (9) Evaluation of the rotational stability of the SISA (Cutting Edge) monofocal intraocular lens after cataract surgery - a prospective uncontrolled open-label study. ID-RCB number: 2017-AOO924-49. Data on file.