BAUSCH + LOMB

LUXSMART™

TORIC

- BG Предварително заредена хидрофобна асферична торична леща
- [Инструкции за потребител \(BG\) – LuxSmart™ Toric](#)
- CS Hydrofobní torické čočky s odděleným zavaděčem - [Navod k pouziti \(CS\) – LuxSmart™ Toric](#)
- DA Forhåndsindsat hydrofob asfærisk torisk linse - [Brugerinstruktioner \(DA\) - LuxSmart™ Toric](#)
- DE Vorgeladene hydrophobe asphärische torische Linse
- [Gebrauchsanweisung \(DE\) – LuxSmart LuxSmart™ Toric](#)
- EL Προφορτωμένος υδροφобος φακός торικός ασφαιρικός - [Οδηγισο χρηστη \(EL\) – LuxSmart™ Toric](#)
- EN Aspheric Toric preloaded hydrophobic lens – [Instructions for use \(EN\) – LuxSmart™ Toric](#)
- ES Lente tórica asférica hidrofóbica precargada - [Instrucciones de uso \(ES\) – LuxSmart™ Toric](#)
- ET Eellaaditud hüdrofoobne asfääriline tooriline lääts - [Kasutusjuhised \(ET\) - LuxSmart™ Toric](#)
- FI Valmiiksi asetettu hydrofobinen asfäärinen toorinen linssi - [Käyttöohje \(FI\) - LuxSmart™ Toric](#)
- FR Lentille hydrophobe asphérique torique préchargée – [Notice d'utilisation \(FR\) – LuxSmart™ Toric](#)
- HR Tvornička hidrofobična asferična tvornička leća - [Upute za upotrebu \(HR\) – LuxSmart™ Toric](#)
- HU Előretöltött hidrofób aszférikus tórikus lencsék - [Hasznalati utasítás \(HU\) – LuxSmart™ Toric](#)
- IT Lente torica asferica idrofobica precaricata - [Istruzioni per l'uso \(IT\) – LuxSmart™ Toric](#)
- LT Asferinis torinis iš anksto įdėtas hidrofobinis lęšis - [Naudojimo instrukcijos \(LT\) - LuxSmart™ Toric](#)
- LV Saspriegta hidrofoba asfēriska toriskā lēca - [Lietotāja instrukcijas \(LV\) - LuxSmart™ Toric](#)
- NL Voorgeladen hydrofobe asferische torische lens - [Instructies voor gebruik \(NL\) – LuxSmart™ Toric](#)
- NO Asfærisk torisk ferdigladet hydrofob linse - [Brukerinstruksjoner \(NO\) - LuxSmart™ Toric](#)
- PL Soczewka hydrofobowa asferyczna toryczna załadowana fabrycznie
- [Instrukcje użytkowania \(PL\) – LuxSmart™ Toric](#)
- PT Lente tórica asférica hidrofóbica pré-carregada - [Instruções de utilização \(PT\) – LuxSmart™ Toric](#)
- RO Lentilă torică asferică hidrofobă preîncărcată - [Instrucțiuni de utilizare \(RO\) - LuxSmart™ Toric](#)
- RU Асферическая торическая гидрофобная линза предварительно загруженная
- [Инструкция по применению \(RU\) – LuxSmart™ Toric](#)
- SK Predsadená hydrofóbna asférická torická šošovka - [Užívateľské pokyny \(SK\) - LuxSmart™ Toric](#)
- SL Prednastavljena hidrofobna asferična torična leča - [Navodila za uporabo \(SL\) - LuxSmart™ Toric](#)
- SV Förladdad hydrofob asfærisk torisk lins - [Användarinstruktioner \(SV\) - LuxSmart™ Toric](#)
- TR Asferik torik önceden yüklenmiş hidrofobik lens - [Kullanici talimatları \(TR\) - LuxSmart™ Toric](#)

Cutting Edge S.A.S.
770 rue Alfred Nobel
Immeuble Le Nobel
34000 Montpellier
FRANCE

CE 0197

INSTRUCTION FOR USE

Date of the last revision: 2024/11/04

Aspec-e19651-04

1- Product to which these instructions for use apply


These instructions for use apply to the following products:

Brand Name	Models	Performance characteristics	Basic UDI-DI
LUXSMART TORIC	SMARTTD (LuxSmart™ Crystal TORIC)	Aspheric Toric preloaded hydrophobic lens	3665493LUXGS8K
	YSMARTTD (LuxSmart™ TORIC)		

2- Description

The LuxSmart™ TORIC intraocular lenses (IOLs) are single-piece ultraviolet (UV) absorbing posterior chamber IOLs developed to replace the natural crystalline lens in aphakic adult patients with pre-existing corneal astigmatism. LuxSmart™ TORIC IOLs have aspheric and toric optics. The LuxSmart™ TORIC lenses incorporate a pure refractive optics (PRO) technology, across its entire optical diameter with no diffractive rings. The lens is provided sterile and is designed to be folded prior to insertion in the eye and implantation in the capsular bag. Preloaded LuxSmart™ TORIC IOLs are supplied with a preloaded Delivery System to provide a touch-free method of delivery of the lens into the capsular bag. Preloaded LuxSmart™ TORIC IOLs are supplied in 2 different models – Models SMARTTD or YSMARTTD. See label on the carton box for type of the lens.

3- Physical Characteristics

Lens / Haptic Material	Model SMARTTD: Foldable hydrophobic acrylic copolymer containing a UV absorber. Model YSMARTTD: Foldable hydrophobic acrylic copolymer containing a UV absorber and a blue light filtering chromophore.
Material Characteristics	LuxSmart™ TORIC hydrophobic acrylic IOLs are made from highly purified acrylate/methacrylate/methylstyrene copolymer (99.9%) including a benzotriazole UV absorber (0.1%). Azo dye (<0.01%) is added to the material of Model YSMARTTD. Index of Refraction when wet at 21°C: 1.544 Index of Refraction when in the eye at 35°C: 1.540
Spectral Transmittance	For +20.00 diopter IOL: Model SMARTTD: 10% transmittance at 393.5 nanometers. Model YSMARTTD: 10% transmittance at 396.7 nanometers.  <p style="text-align: center;">with chart's X-value = Wavelength (nm) and Y-value = % Transmittance</p>
Optic Type	Models LuxSmart™ Crystal TORIC and LuxSmart™ TORIC: Aspheric Toric
Powers (SE)	+06.00 to +34.00 Diopters
Cylindrical Powers	+0.75 to +6.00 Diopters
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 11 mm; Haptic Angle: 0°
Lens Orientation	The IOL is to be implanted with the anterior side of the lens facing up towards the anterior side of the eye. The orientation of the IOL can be verified by visual inspection of the lens. As illustrated in FIG.4, when the orientation features are top right (A) and bottom left (B), you are facing the anterior side of the lens.

4- Indications

Intended purpose:

LuxSmart™ TORIC IOLs are intended to replace the human crystalline lens.

Indication:

The LuxSmart™ TORIC posterior chamber lens is indicated for primary implantation for the visual correction of aphakia and pre-existing corneal astigmatism in patients in whom the cataractous lens has been removed. The lens is designed for implantation in the capsular bag.

Patient target group:

Adult patients requiring visual correction of aphakia, regardless of ethnicity or gender with pre-existing regular corneal astigmatism.

Intended users:

IOLs must be handled by health professionals and implanted as part of a surgical procedure by a qualified ophthalmic surgeon.

Mode of action:

LuxSmart™ TORIC IOLs positioned in the capsular bag act as a refractive medium to correct aphakia and pre-existing astigmatism. LuxSmart™ TORIC IOLs have a posterior toric optical surface with axis marks that identify the flat meridian of the cylinder (plus cylinder axis) (FIG.4). Alignment of the axis marks with postoperative steep corneal meridian (intended axis of placement) ensures correct placement and optimal results.

Lifetime:

The LuxSmart™ TORIC IOL is expected to be stable indefinitely over the lifetime of the patient.

5- Clinical benefits

Provide a functional corrected distance visual acuity (≤ 0.3 logMAR) and predictable refractive outcomes.

6- Safety and clinical performance

The summary of safety and clinical performance (SSCP) is published in the European Database on Medical Devices (Eudamed).

The URL of the public website of Eudamed is: <https://ec.europa.eu/tools/eudamed>.

The SSCP is linked to the Basic UDI-DI listed in section 1.

7- Packaging/Sterilization

LuxSmart™ TORIC IOLs are individually packaged preloaded in a delivery system in a blister. All models are supplied dry and should be opened under sterile conditions. A patient card and self-adhesive labels are supplied in the carton to provide traceability of the IOL. LuxSmart™ TORIC IOLs are labeled with spherical equivalent and cylindrical powers. Spherical equivalent corresponds to the power obtained with most calculation formulas.

All models are sterilized by Ethylene Oxide.

8- Precautions for use and storage

- Store at room temperature. Avoid temperatures below 5°C and above 45°C. Do not freeze.
- LuxSmart™ TORIC IOLs, balanced salt solution, and the viscoelastic solution used should be stored for at least 30 min at minimum 20°C before surgery.
- A high level of surgical skill is required to implant intraocular lenses. Surgeon should be fully trained before attempting to implant intraocular lenses.
- The IOL model, its powers and expiration date should be verified before opening the protective packaging and before opening the individual sterile blister.
- Sterility of the device is guaranteed only if the individual sterile blister has not been opened or damaged.
- Handle the lens carefully to avoid surface or haptics damages.
- The IOL must be carefully rinsed prior to implantation. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution.
- Care should be taken to verify that the lens is fully implanted into the capsular bag.
- If a YAG laser posterior capsulotomy is performed, assure that the laser beam is focused slightly behind the posterior capsule.
- Customization of surgically induced astigmatism (SIA) is recommended based upon individual surgical technique and past results since size and location of the surgical incision may affect amount and axis of postoperative corneal astigmatism.
- Consideration of posterior corneal astigmatism is recommended for optimal results.
- Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens after implantation in the capsular bag, not to overinflate the capsular bag and to ensure that the wound is watertight at the end of the case. Residual viscoelastic / overinflation / wound leak may allow the lens to rotate away from the intended axis of placement.
- Special care should be taken to ensure proper orientation of the toric lens at the intended axis of placement at the end of the surgical case.

- Consideration of pupil size is recommended as sufficiently dilatable pupils are of importance to visualize axis marks.

Disposal

The injection system and discarded IOLs (used or unused) are classified as medical (clinical) waste that harbors a potential infection or microbial hazard and must be disposed of accordingly.

Disposal should be done in compliance with the local regulation of waste from health-care activities characterized by risk of infection.

9- Directions for use

1. The physician should determine preoperatively the powers of the lens to be implanted, using calculation methods and accurate biometry measurements depending on the surgeon's experience, preferences and/or with the assistance of a Toric Calculator. The recommended A-constants listed on the lens carton are guidelines only.

2. Prior to surgery, using the physician's preferred method, the operated eye should be referenced to provide a clear indication of the intended axis of placement.

3. Prior to implanting, examine the lens package for model, powers (spherical equivalent and cylinder), expiration date and proper configuration.

4. Open the carton and inspect the sterilized blister containing the lens and the injection system. Make sure it is not damaged. Verify that the information on the blister label is consistent with the information on the carton box and the self-adhesive labels. Peel open the blister and transfer the device onto the sterile field.

5. Inspect the device and make sure it is intact. The lens should be rinsed thoroughly by flushing balanced salt solution through the tip of the cartridge (FIG. 1a) and from the end of the loading chamber (FIG. 1b). Full hydration of the cartridge and the loading chamber is mandatory.

Please refer to the enclosed figures to fully prepare the device before implantation, taking care not to touch the IOL:

FIG. 1a: Apply viscoelastic only in the tip of the cartridge. Wait for 30 seconds to ensure lubrication (coating activation).

FIG. 2: Close the wings of the loading chamber. When the "click-lock" mechanism is engaged, the lens is securely in place. Preparation is complete.

6. Press the injector plunger forward until the plunger stops in the middle of the viewing window (FIG. 3).

The lens should be positioned in the viewing window until implantation and should be injected in the shortest possible time after closing the loading chamber.

7. Insert the cartridge tip in the incision with the bevel of the injector tip facing down to the patient's eye.


In push mode: open the finger flanges and inject by pushing the plunger forward smoothly but continuously.

In screw mode: inject by rotating the screw mechanism clockwise to advance the lens forward until the lens is delivered into the capsular bag.

8. Discard the injection device (See §8 Disposal).

9. The lens is to be implanted with the anterior side of the lens facing up towards the anterior side of the eye. The orientation of the IOL can be verified by visual inspection of the haptics. As illustrated in FIG. 4, when the haptic insertions are top right (A) and bottom left (B), you are looking at the anterior side of the lens.

10. Align LuxSmart™ TORIC IOL's axis marks with the referenced intended axis of placement.

There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient.

10- Contra-indications and Warnings

Physicians considering lens implantation in any of the following circumstances (but not limited to) must weigh the potential risk/ benefit ratio: recurrent severe anterior or posterior segment inflammation or uveitis; patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases; surgical difficulties increasing the potential for complications (such as but not limited to: persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss); distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible; circumstances that would result in damage to the endothelium during implantation; patients in whom neither the posterior capsule nor zonule are intact enough to provide support; suspected microbial infection. Toric IOLs are not recommended with pre-existing irregular astigmatism, radial tears on the capsulotomy, zonular damage or when the need for a large capsulotomy is anticipated.

Warnings

The IOL must be implanted in accordance with the instructions for use. Improper use may impose risk to the patient's health.

Do not use if product sterility or quality is thought to be compromised due to damaged packaging or if product was unintentionally opened before use.

Do not attempt to modify or alter the device or any of its components.

Do not use after expiration date.

Do not use without viscoelastic solution.

Do not reuse. Reuse may impose serious risk to the patient's health.

Do not resterilize as this can produce undesirable side effects.

Do not implant the lens in the anterior chamber or in the ciliary sulcus.

Rotation of LuxSmart™ TORIC lens away from its intended axis can reduce its astigmatic correction.

Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, repositioning should occur as early as possible prior to lens encapsulation.

Do not use refractive data for calculation since the presence of lenticular astigmatism in the natural lens to be removed may influence results.

Amvisc® Plus, OcuCoat® or EyeFill® C viscoelastic should be used.

Only the injection system provided with the preloaded LuxSmart™ TORIC lens should be used.

11- Possible undesirable side effect

As with any surgical procedure, there is risk involved. Potential complications of cataract surgery may include, but are not limited to the following: secondary cataract, inflammation (iridocyclitis, pupillary/cyclitic membrane, vitreous inflammation, CME, TASS ...), elevated IOP, anterior capsule opacification (fibrosis, phimosis), corneal endothelial damage, endophthalmitis, retinal detachment, corneal edema, pupillary block, iris prolapse, hypopyon, glaucoma, wound leak, hyphema, iris atrophy, vitreous prolapse, IOL decentration or tilt, IOL optic deposits or opacities, glistenings, residual refractive error, visual symptoms and secondary surgical intervention (such as but not limited to wound leak repair, vitreous aspiration, iridectomy, IOL repositioning or replacement and retinal detachment repair).

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this IOL or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

12- Implant card and patient information




The implant card included in the package is to be completed and given to the patient together with instructions to keep this card as a permanent record of the implant and to show the card to any eye care professional consulted in the future.

Patient information is made available on the internet. The link to access the information is printed on the implant card.

13- Information to be conveyed to the patient

- Patients should be advised that unexpected outcomes could lead to continued spectacle dependence.

14- Symbol used on labelling

	Manufacturer		Do not re-use
	Date of manufacture		Consult instructions for use or consult electronic instructions for use
	Use-by date		Caution
	Batch code		Patient name
	Catalogue number		Health care center or doctor
	Serial number		Date
	Sterilized using ethylene oxide		Medical device
	Do not re-sterilize		Unique device identifier
	Do not use if package is damaged and consult instructions for use	\emptyset_B	Diameter of the body
	Single sterile barrier system with protective packaging outside	\emptyset_T	Overall diameter
	Lower limit of temperature	D	Diopter
	Right eye		Left eye

 Cutting Edge S.A.S.
770 rue Alfred Nobel
Immeuble Le Nobel
34000 Montpellier
FRANCE

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